RaySafe X2 Solo DENT

X-ray meter for dental applications





Quality Assurance in Dental X-ray


Hundreds of millions of X-ray exposures are done every year in dental clinics around the globe. The majority are intraoral where the dose is relatively low, but these machines need to be checked, following the ALARA principles. Globally there are lots of analogue intraoral machines that may give many times the necessary dose, which is easily checked using a meter from RaySafe.

At your fingertips

From a measurement perspective, intraoral machines are easy to check. Panoramic, however is more challenging due to a very narrow radiation field, which requires accurate positioning of the sensor. RaySafe X2 Solo DENT makes these measurements easier than ever and is your complete tool for all measurement needs in the dental area.

Simplicity has always been a hallmark for Unfors RaySafe. Our simplicity comes from a careful balance of advanced sensor technology and intelligent signal processing.

The result – a device that provides accurate measurements with the ultimate in user friend-liness. Right at your fingertips.



RaySafe X2 Solo DENT

- Leading precision designed for specific needs

The RaySafe X2 Solo DENT is a meter for quality assurance and service of dental X-ray machines. This fast and easy-to-use meter provides leading precision and tailor made measurement capabilities for all dental applications. RaySafe X2 Solo DENT covers all your needs:

- Supports all types of dental machines;
 Cone Beam CT, Panoramic and Intraoral.
- Options include HVL & Total Filtration as well as invasive mAs measurements.
- Measures dose, dose, rate, kVp, exposure time, pulses, pulse rate and dose/pulse in the kV and dose rate ranges needed for dental applications.

In addition, with its minimal setup time, acquiring your first exposures take less than one minute. Because of its intuitive user interface and built-in intelligence, you can rapidly learn how to use the RaySafe X2 Solo DENT. The user will be able to focus on interpreting the measured data instead of focusing on how to obtain data. Ease-of-use saves time, but most importantly, it minimizes the risk of costly user errors.

Panoramic measurements

RaySafe X2 Solo DENT includes the brand new X2 Panoramic Holder. Correctly positioning the sensor in the narrow X-ray field of a Panoramic machine is now easier than ever:

- 1. Attach the holder to the X-ray machine.
- 2. Make an exposure with gafchromic film (included) to reveal the position of the X-ray beam.
- 3. Use the adjustment lever to align the sensor with the X-ray beam.

The markings on the holder makes it possible to go back to a previous position if needed.



X2 PANORAMIC HOLDER

RaySafe X2 technology

RaySafe X2 Solo DENT is built on the RaySafe X2 platform, highly esteemed for its user-friendliness and performance. Features of the X2 Solo DENT include:



WAVEFORM Overview and simple analysis of kVp, dose rate

LARGE TOUCH SCREEN

display for simple operation and a great overview of all measured parameters.

FULL WAVEFORMS

directly in the base unit for quick analysis of measurements.

NO SPECIAL SETTINGS

to handle different types of X-ray machines. Just connect and measure.

STACKED SENSORS

for easy positioning and minimum radiological footprint.

WIDE DYNAMIC RANGE

 no special modes for high sensitivity measurements.

BUILT-IN MEMORY

 up to 10 000 measurements with waveforms are stored in the base unit.

Options for RaySafe X2 Solo DENT

The following options are available to the RaySafe X2 Solo DENT:

HVI & TOTAL FILTRATION

With this option, RaySafe X2 Solo DENT will measure Half Value Layer (HVL) and Total Filtration in a single exposure avoiding the need of multiple dose measurements with different filters. RaySafe X2 Solo utilizes a stacked sensor with diodes having different filtration, and hence accurately calculates HVL. With the HVL & Total Filtration option installed, these parameters will be measured simultaneously as the other parameters and will be displayed on the main screen after each exposure.

RAYSAFE SERVICE PROGRAM

The RaySafe Service Program makes it possible to control your life cycle cost and will keep your instrument in perfect functional condition year after year. Available with a 12 or 24 months cycle, the RaySafe service program will extend the warranty of your meter up to 8 years. Also, members of the service program are prioritized and hence enjoy a faster turnaround time. The RaySafe service program can be pre-paid at the time of instrument purchase.

mAs

RaySafe X2 Solo DENT is available with or without mAs. The design allows to invasively measure tube current even for pulsed measurements.

X2 SOLO GENERAL

EMCAccording to IEC 61326-1SAFETYAccording to IEC 61010-1X-RAY METERS STANDARDComplies with IEC 61674

EXPOSURES NEEDED One

USB CABLES 2 m (6.6 ft), 5 m (16.4 ft) and

5 m active extender

SIZE BASE UNIT 34 x 85 x 154 mm (1.3 x 3.3 x 6.1 in)

WEIGHT BASE UNIT 521 g (18.4 oz)

OPERATING TEMPERATURE 15 – 35 °C (59 – 95 °F)

STORAGE TEMPERATURE -25 – 70 °C (-13 – 158 °F)

POWER SOURCE Rechargeable Li ion battery

BATTERY TIME ~10 hours intensive usage

BATTERY TESTED According to UN 38.3

DISPLAY 4.3" LCD with capacitive touch

MEMORY ~10000 latest exposures


SOFTWARE X2 View for data handling and

analysis. Also exports data to

Microsoft Excel.

TAILOR MADE FOR DENTAL MEASUREMENTS

The RaySafe X2 Solo DENT provides high precision measurements of all essential X-ray parameters for dental applications; Cone Beam CT, Panoramic and Intraoral X-ray. RaySafe X2 Solo DENT is designed to improve your productivity by its intuitive user interface, built-in intelligence, and minimal set-up time.



X2 mAs

RESOLUTION

mAs

RANGE 0.001 – 9999 mAs

0.001 mAs

UNCERTAINTY 1%

mΑ

RANGE (PEAK) 0.1 – 1500 mA

RESOLUTION 0.01 mA

UNCERTAINTY 1%

TIME

RANGE 1 ms – 999 s

RESOLUTION 0.1 ms **BANDWIDTH** 1 kHz

UNCERTAINTY 0.5 %

PULSES

RANGE 1 – 9999 pulses

RESOLUTION 1 pulse

PULSE RATE

RANGE 0.1 – 200 pulses/s

RESOLUTION 0.1 pulse/s

mAs/PULSE

RANGE 0.001 – 9999 mAs

RESOLUTION 0.001 mAs

UNCERTAINTY 1%

WAVEFORM

RESOLUTION 125 μ s* **BANDWIDTH** 1 kHz

^{*} automatically reduced for exposures longer than 3 s

X2 DENT SENSOR

WEIGHT 42 g (1.5 oz)

SIZE 14 x 22 x 79 mm (0.5 x 0.9 x 3.1 in)

ACTIVE COMPENSATION

Beam quality independent for the following ranges:

 DOSE/DOSE RATE
 40 – 130 kVp, 1 – 14 mm AI HVL

 kVp
 40 – 130 kVp, up to 1 mm Cu

 TF
 60 – 120 kVp, up to 1 mm Cu

DOSE

RANGE 1 nGy – 9999 Gy (0.1 μR – 9999 R)

UNCERTAINTY 5 % or 5 nGy (0.5 μ R)

DOSE RATE

 $\textbf{RANGE} \hspace{1.5cm} 1\,\mu\text{Gy/s} - 500 \text{ mGy/s}$

(5 mR/min – 3400 R/min)

 RESOLUTION
 1 nGy/s (5 μR/min)

 TRIG LEVEL
 1 μGy/s (7 mR/min)

UNCERTAINTY 5 %

kVp

RANGE 40 - 130 kVp **MINIMUM DOSE** $50 \mu\text{Gy (6 mR)}$

MINIMUM DOSE RATE (PEAK) 10 μGy/s (70 mR/min)

UNCERTAINTY 2 %

HVL (OPTIONAL)

 RANGE
 1 – 14 mm Al

 MINIMUM DOSE
 1 μGy (120 μR)

MINIMUM DOSE RATE (PEAK) $1 \mu Gy/s$ (7 mR/min) at > 70 kV

 $2.5~\mu\text{Gy/s}$ (17 mR/min) at 50 kV

UNCERTAINTY 10 %

TOTAL FILTRATION (OPTIONAL)

 RANGE
 1.5 – 35 mm Al

 MINIMUM DOSE
 50 μGy (6 mR)

 MINIMUM DOSE RATE (PEAK)
 10 μGy/s (70 mR/min)

 UNCERTAINTY
 10 % or 0.3 mm Al

TIME

RESOLUTION

RANGE 1 ms - 999 s

BANDWIDTH 4 Hz – 4 kHz*

0.1 ms

UNCERTAINTY 0.5 %

* automatically adjusted depending on signal level

PULSES

RANGE 1 – 9999 pulses

PULSE RATE

RANGE 0.1 – 200 pulses/s

DOSE/PULSE

RANGE 1 nGy/pulse – 999 Gy/pulse

(0.1 μR/pulse – 999 R/pulse)

WAVEFORMS

RESOLUTION 62.5 μ s*

BANDWIDTH kV $0.1 - 0.4 \text{ kHz}^{**}$

BANDWIDTH DOSE RATE 4 Hz – 4 kHz**

UNFORS RAYSAFE UNCERTAINTY DEFINITION

The expanded uncertainty is stated as the combined uncertainty of measurement multiplied by the coverage factor k=2, which assuming a normal distribution has a coverage probability of 95 % (complies with GUM by ISO (1995, ISBN 92-67-10188-9)).

Instrument specifications are subject to purchased configuration. All specifications may change without notice.

^{*} automatically reduced for exposures longer than 1.5 s

^{**} automatically adjusted depending on signal level

ADDITIONAL DENTAL X-RAY QA SOLUTIONS







RAYSAFE THINX INTRA

RAYSAFE X2

PRO-DENT SET

Unfors RaySafe offers comprehensive solutions for the X-ray room to measure the performance of X-ray equipment and to monitor medical staff dose in real-time.

RaySafe helps you avoid unnecessary radiation.

